

# IMAGE-COM<sup>©</sup>

powered by  **IMAGE ARENA**



**CARDIOLOGY**

Comprehensive Image-Review,  
Measurements and Reporting





# IMAGE-COM<sup>©</sup>


## Comprehensive Image-Review, Measurements and Reporting

Time efficient Reporting and Image-Review plays a decisive role in improving the quality of patient care. Image-Com<sup>©</sup> has been developed to speed up workflow in the daily routine. It makes cardiological information available faster, provides additional clinical information and accelerates reporting.

### Image-Com<sup>©</sup>\*

As a viewer for cardiac catheter and cardiac ultrasound examinations Image-Com<sup>©</sup> offers a comprehensive feature set for rapid and clearly structured review of still images and image sequences.


- Fast and flexible composition of relevant image data ensured by intuitive user interface
- Side-by-side comparison of current and previous examinations (XA and Echo)
- Comfortable analysis of image sequences using video functions
- Export of image-data to AVI, BMP, JPEG or DICOM format



### Image-Arena Measurement\*

Image-Arena<sup>™</sup> Measurement provides measurement and re-measurement capabilities for Echo examinations.


- Automatic import of measurement values from ultrasound systems
- Availability of all current 2D Echo-, M-Mode- and Doppler measurements
- Measurement of up to five values for each parameter



### Image-Arena Reporting\*

Based on report templates Image-Arena<sup>™</sup> Reporting ensures accelerated generation of standardized Echo reports

- Automatic import of measurement values
- Rapid report text generation using text modules
- Integration of diagnosis relevant images by drag & drop



\* Optional module of Image-Arena<sup>™</sup>  
Image-Arena<sup>™</sup> is a trademark of TomTec Imaging Systems GmbH

