



SONOLiver[®]

TOMTEC
IMAGING SYSTEMS




Benign or Malignant?



Increase Diagnostic Confidence



For a higher specificity in tumor characterization, *sonoLiver*® provides unique features for the perfusion quantification of focal liver lesions (FLL) by using dynamic vascular pattern (DVP) display. *sonoLiver*® offers objective parameters to increase diagnostic confidence.



Clinical Use

- Improved characterization of suspicious lesions
- Better differentiation between benign and malignant lesion types
- Increased diagnostic confidence

Major Benefits

- Convenient documentation by automatic report generation
- Improved diagnostic reliability by enhanced visualization of DVPs in comparison to healthy parenchyma
- System independent results by linearization of data of different ultra sound devices



Images by courtesy of J.-Y. Meuwly MD, University Hospital Lausanne.

sonoLiver® is only available for sale in specific countries. For availability, please contact your local representative. *sonoLiver*® is a software distributed by TomTec Imaging Systems GmbH and a registered trademark of Bracco International BV.

